The Ahmed[®] Glaucoma Valve

Model FP7

NEW WORLD
MEDICAL

CE
0459

Index

English (en)	3
Français (fr)	9
Deutsch (de)	15
Italiano (it)	21
Español (es)	27
Português (pt)	33
Čeština (cs)	39
Dansk (da)	45
Ελληνικά (el)	51
Suomi (fi)	57
Hrvatski (hr)	63
Latviski (lv)	69
Magyar (hu)	75
Nederlands (nl)	81
Polski (pl)	87
Română (ro)	93
Slovenščina (sl)	99
Svenska (sv)	105
Lietuvių (lt)	111
Norsk (no)	117
Íslenska (is)	123

The Ahmed® Glaucoma Valve Model FP7

Instructions for Use

DESCRIPTION:

The Ahmed® Glaucoma Valve Model FP7 (AGV-FP7) is a valved aqueous drainage implant designed to regulate intraocular pressure in eyes suffering from refractory glaucoma. The Ahmed® device consists of a valve mechanism and a plate. The valve mechanism is comprised of a silicone drainage tube and polypropylene body which houses a silicone elastomer valve membrane. The polypropylene body protects the valve membrane from blockage by fibrous tissue. The plate conforms to the shape of the globe at its equator and provides a surface where fluid can be dispersed. The plate material of the Model FP7 is silicone. The device is for single use only. Once implanted, the Ahmed® Glaucoma Valve remains permanently in the patient's eye. The device is expected to be stable over the lifetime of the patient.

INDICATIONS:

The Ahmed® Glaucoma Valve Model FP7 is indicated for the management of refractory glaucomas, where previous surgical treatment has failed, or by experience is known not to provide satisfactory results. Such refractory glaucomas can include neovascular glaucoma, primary open angle glaucoma unresponsive to medication, congenital or infantile glaucoma, and refractory glaucomas resulting from aphakia or uveitis.

INTENDED PATIENT POPULATION:

• Adult and pediatric patients.

CLINICAL BENEFITS:

Immediate reduction in intraocular pressure (IOP) and/or reduction in the need for IOP lowering medications. IOP is the only modifiable risk factor for glaucoma progression. Glaucoma progression can lead to blindness. These clinical benefits may diminish over time due to the natural wound healing response around the implant which restricts fluid outflow to the surrounding tissue. The valve mechanism reduces the likelihood of hypotony compared to non-valved glaucoma drainage devices which do not have a built-in outflow restriction mechanism.

Consistent with the European Glaucoma Society (EGS) guidelines: "there is no single target IOP level that is appropriate for every patient, so the target

IOP needs to be estimated separately for each eye of every patient and adapted to changing patient status". Further, the EGS states "Target IOP should be updated at each monitoring visit on the basis of changes in glaucoma or other ocular or systemic diseases."

The document provides factors to consider when setting the target IOP. "There is little evidence to support any particular algorithm to set the target IOP. In newly diagnosed patients, the target IOP is initially determined according to the stage of disease and the baseline IOP. The treatment goal is typically estimated as a specific pressure level or a percentage reduction." It is up to the healthcare professional to judge the level of IOP sufficient to prevent the risk of progression of sight loss.

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE:

The summary of safety and clinical performance (SSCP) for Ahmed Glaucoma Valve is available in the European database on medical devices (EUDAMED) and it is linked to the basic UDI-DI: 00892064002119.

<https://ec.europa.eu/tools/eudamed>

CONTRAINDICATIONS:

- | | |
|-----------------------------|-----------------------------|
| a. Bacterial conjunctivitis | e. Bacteremia or septicemia |
| b. Bacterial corneal ulcers | f. Active scleritis |
| c. Endophthalmitis | g. No light perception |
| d. Orbital Cellulitis | |

COMPLICATIONS AND ADVERSE REACTIONS:

Complications and adverse reactions during or following surgery may include:

- | | |
|----------------------|----------------------|
| Corneal Edema | Choroidal detachment |
| Corneal Touch | Iritis |
| Iris/Tube Touch | Hyphema |
| Synechia | Tube Obstruction |
| Exposed Sclera Graft | Tube Retraction |

as well as known complications of aqueous shunts and general intraocular surgery including hypotony, shallow chamber, hyphema, choroidal effusion, suprachoroidal hemorrhage, choroidal detachment, retinal detachment, iritis, synechia, cataract, conjunctival buttonhole, phthisis bulbi, bullous keratopathy, uveitis and endophthalmitis.

WARNINGS, PRECAUTIONS:

Do not use the device if sterile package integrity has been compromised. Do not re-sterilize the implant. Before using the Ahmed[®] Glaucoma Valve, the implanting surgeon should be skilled in glaucoma filtering procedures and familiar with the use of drainage devices, as well as post-operative care required.

Priming the valve is essential for proper functioning of the valve. Do not implant the valve without priming it.

Ensure that the valve is primed before operating on the patient. Inherently, the force required to prime the valve may vary between units. Some valves may require more force applied over a longer period of time than other valves in order to perform the priming operation. If you are unable to prime the valve, do not use the device. The manufacturer will exchange the device.

In rare cases, you may be unable to prime the valve due to crimped or pinched tubing. In such a situation, do not use the device. The manufacturer will exchange the device.

Using a needle for priming can puncture the tube resulting in an undesirable leak or inability to prime. Only use a blunt cannula to prime the valve.

The Ahmed[®] Glaucoma Valve is intended for single use, it shall not be reused. It shall be discarded to avoid deterioration, cross contamination and/or infection.

Tampering with the valve can cause malfunctioning of the valve.

STORAGE CONDITIONS:

Store Ahmed[®] Glaucoma Valve at room temperature.

Caution: Do not expose to direct sunlight and/or water.

PATIENT INFORMATION:

The surgeon should give the patient the implant card provided with the appropriate information filled in and should advise the patient to keep the card in a safe place, e.g., his or her wallet, for future reference. The surgeon should advise the patient that this implant card contains important information related to the Ahmed[®] Glaucoma Valve implant and that the card should be

shown to their current and future health care providers. The surgeon should inform patient the implant card contains the following information:

- Implant date
- Name and Address of the implanting healthcare institution / provider
- Serial number
- Lot number
- The Unique Device Identifier (UDI)
- Model number
- Address and website of the manufacturer
- Website for patient information

The surgeon should advise the patient that the Ahmed[®] Glaucoma Valve is comprised of materials with a long history of safe clinical use that remain stable within the ocular environment over the lifetime of the patient.

The surgeon should provide information to the patient that there is no device specific post-op care outside of normal generic ophthalmic instructions.


REPORTING:

Adverse reactions and/or potentially sight-threatening complications that may be reasonably regarded as product related and that were not previously expected in nature, severity or incidence must be reported to New World Medical (NWM). Report any U.S. product related event to the toll free NWM number at +1 (800) 832-5327. For customers from countries outside of the U.S. please contact your local company representative. For customers from countries within the EU, any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

AGV-FP7

Trade Name: Ahmed® Glaucoma Valve Model FP7

Top View

Side View

End Plate: Medical Grade Silicone
Drainage Tube: Medical Grade Silicone
Valve Membrane: Medical Grade Silicone
Valve Casing: Medical Grade Polypropylene (PP)
Adhesive: Medical Grade Silicone Glue

Maximum Thickness: 2.1 mm
Max Length: 41.4 mm
Plate Length: 16.0 mm
Device Overall Width: 13.0 mm
Drainage Area/Explant Surface Area: 184 mm²

The implant is supplied sterile in a sealed pouch. Product Information, Patient Chart Labels, Implant Card, and Implant Card Instructions are also enclosed with the sterile package. The implant has been terminally sterilized by Gamma Radiation. Sterility is assured provided that the peel pouch has not been opened or damaged, and the sterility expiration date has not lapsed. The peel pouch is intended to be opened using sterile technique, allowing the implant to be dropped into the sterile field.

OPENING DEVICE AND MAINTAINING STERILE FIELD:

To open packaging and perform device transfer to the surgeon and sterile field, the following sequence of steps should be performed.

- Verify the model number to be used and that the expiration date has not passed.
- Open the carton, remove, and inspect the contents for damage using aseptic techniques.
- Remove the inner pouch by peeling open the outer pouch using the flap adjacent to the chevron shaped seal using aseptic techniques.
- Repeat the process for opening the inner pouch and remove the implant using aseptic techniques.

- The surgeon should examine the implant prior to use for integrity and to ensure particles are not present, as it may become electrostatically charged upon opening of the package. Rinse the device in sterile saline if required.

DISPOSAL:

1. Dispose of explanted devices in properly labeled Biohazard waste containers.
2. The Ahmed® Glaucoma Valve is a single patient, single use device. Do not reuse or re-sterilize. Re-sterilized or reused instruments pose risks of infection and may result in damage to tissue.

SURGICAL PROCEDURE

Figure 1

Step 1


Step 2


Step 3


Step 4


Step 5

Step 6

Step 7

Alternative to Step 7

Refer to figure 1 for the illustration of each step. The steps described here are intended as a guideline only, and do not represent recommended treatment for any particular patient. The use of any specific surgical technique or maneuver is at the sole discretion of the surgeon. Surgeons should be familiar with the use of glaucoma drainage devices and post-operative care considerations before implanting any drainage device. Reference papers and surgical video tapes are available upon request.

1. The implant should be examined and primed prior to implantation. Priming is accomplished by injecting 1 cc balanced salt solution or sterile water through the drainage tube and valve, using a blunt 26-gauge cannula.
2. A fornix-based incision is made through the conjunctiva and Tenon's capsule. A pocket is formed at the superior quadrant between the medial or lateral rectus muscles by blunt dissection of Tenon's capsule from the episclera.
3. The valve body is inserted into the pocket between the rectus muscles and sutured to the episclera. The leading edge of the plate should be at least 8-10 mm from the limbus.

4. The drainage tube is trimmed to permit a 2-3 mm insertion of the tube into the anterior chamber (AC). The tube should be bevel cut to an anterior angle of 30° to facilitate insertion.
5. A paracentesis is performed, and the AC is entered at 1-2 mm way from the limbus with a sharp 23-gauge needle to create a needle track, parallel to the iris. Caution: Care must be taken to ensure that the drainage tube does not contact the iris or corneal endothelium after insertion.

Note: Some surgeons prefer to enter the AC from at least 3 mm away from the limbus.

6. The drainage tube is inserted approximately 2-3 mm into the AC through the needle track created in step 5. The leading edge of the plate should be 8-10 mm from the limbus.
7. The exposed drainage tube is covered with a piece of preserved, donor sclera, pericardium, cornea, or other suitable patch graft material which is sutured into place and the conjunctiva is closed.

Note: As an alternative to Step 7, a 2/3 thickness limbal-based scleral flap may be made. The tube is inserted into the AC through a 23-gauge needle puncture made under the flap. The flap is sutured closed.

EXPLANT:

If removal of the device is necessary, standard surgical technique can be used to isolate the plate and tube from surrounding conjunctiva and Tenon's capsule.

1. Any adhesions between the device and surrounding tissues should be amputated.
2. Pre-existing sutures are cut and removed to allow for free mobility of the plate and tube.
3. The silicone tube is then grasped and removed from the sclerotomy tunnel.
4. The entire device is then removed from the operative site.
5. The sclerotomy tunnel is sutured closed and watertight followed by re-approximation of the conjunctiva and Tenon's capsule using the surgeon's preferred suture and technique.

CLINICAL STUDY RESULTS:

Clinical data for the AGV-FP7 can be found in the citations referenced below. This includes safety and effectiveness data and information on adverse effects and complications associated with the AGV-FP7.

1. The Ahmed Baerveldt Comparison Study. This study compared the long-term outcomes and complication of the Ahmed Glaucoma Valve model FP7 (AGV-FP7) and the Baerveldt glaucoma implant. This was a multicenter, prospective, randomized, controlled clinical trial with several relevant publications which have been cited here.

Barton K, Gedde SJ, Budenz DL, Feuer WJ, Schiffman J. The Ahmed Baerveldt Comparison Study methodology, baseline patient characteristics, and intraoperative complications. *Ophthalmology*. 2011 Mar;118(3):435-42. PubMed PMID: 20932581; NIHMSID: NIHMS233027; PubMed Central PMCID: PMC3020244.

Budenz DL, Barton K, Feuer WJ, Schiffman J, Costa VP, et al. Treatment outcomes in the Ahmed Baerveldt Comparison Study after 1 year of follow-up. *Ophthalmology*. 2011 Mar;118(3):443-52. PubMed PMID: 20932583; NIHMSID: NIHMS233028; PubMed Central PMCID: PMC3020266.

Budenz DL, Barton K, Gedde SJ, Feuer WJ, Schiffman J, et al. Five-year treatment outcomes in the Ahmed Baerveldt comparison study. *Ophthalmology*. 2015 Feb;122(2):308-16. PubMed PMID: 25439606; NIHMSID: NIHMS625185; PubMed Central PMCID: PMC4306613.

Budenz DL, Feuer WJ, Barton K, Schiffman J, Costa VP, et al. Postoperative Complications in the Ahmed Baerveldt Comparison Study During Five Years of Follow-up. *Am J Ophthalmol*. 2016 Mar;163:75-82. PubMed PMID: 26596400; NIHMSID: NIHMS746915; PubMed Central PMCID: PMC4901387.

2. Comparison of Polypropylene and Silicone Ahmed Glaucoma Valves. This study was conducted to evaluate and compare clinical outcomes after implantation of the silicone plate Ahmed Glaucoma Valve Model FP7 (AGV-FP7) and the polypropylene plate Ahmed Glaucoma Valve Model S2 (AGV-S2). This was a prospective, multicenter, comparative study.

Ishida K, Netland PA, Costa VP, Shiroma L, Khan B, et al. Comparison of polypropylene and silicone Ahmed Glaucoma Valves. *Ophthalmology*. 2006 Aug;113(8):1320-6. PubMed PMID: 16877071.

3. Comparison of silicone and polypropylene Ahmed glaucoma valves: two-year follow-up. The purpose of this study was to compare the safety and efficacy of the silicone Ahmed glaucoma drainage device (AGV-FP7) and the polypropylene Ahmed glaucoma drainage device (AGV-S2). This study was a single surgeon, retrospective, consecutive case series.

Mackenzie PJ, Schertzer RM, Isbister CM. Comparison of silicone and polypropylene Ahmed glaucoma valves: two-year follow-up. *Can J Ophthalmol*. 2007 Apr;42(2):227-32. PubMed PMID: 17392844.

RETURN OF DAMAGED / UNSTERILE PRODUCT:

Contact NWM for product return policy. Return the device with proper identification and reason for return. Label the return as a biohazard.

Contact Number: US: +1 (909) 466-4304

Email: orders@newworldmedical.com

Return Address:

New World Medical
10763 Edison Court
Rancho Cucamonga
California 91730, USA

Symbols Used on Packaging:

	Manufacturer		Authorized representative in the European Community		Sterilized using irradiation
	Date of manufacture		Consult instructions for use or consult electronic instructions for use		Do not use if package is damaged
	Use-by date		Read all warnings and precautions in instructions for use		MR safe
	Serial number		Do not reuse		Quantity of Device in Package
	Batch Code		Do not re-sterilize		Catalogue Number
	Single Sterile Barrier system with protective packaging outside		Distributor		Medical Device
	Double Sterile Barrier System		Keep Dry		Keep away from sunlight
	Single Sterile Barrier System with Protective packaging inside		Caution: Federal (USA) law restricts this device to sale by or on the order of a physician		Unique Device Identifier
	Open Here				

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California, 91730 U.S.A.
Emergo Europe B.V.