The Ahmed ClearPath® Glaucoma Drainage Device

Model CP250 and CP350

NEW WORLD
MEDICAL

CE
0459

Index

English (en)	3
Español (es).	9
Français (fr)	15
Italiano (it)	21
Nederlands (nl).	27
Português (pt).	33
Dansk (da)	39
Deutsch (de).	45
Ελληνικά (el)	51
Suomi (fi)	57
Polski (pl)	63
Svenska (sv)	69
Norsk (no)	75

The Ahmed ClearPath® Glaucoma Drainage Device

Figure 1


Figure 2


Figure 3


Figure 4


Figure 5



Figure 6


Device Description:

The Ahmed ClearPath® glaucoma drainage device is a non-valved drainage device designed to shunt aqueous in eyes suffering from refractory glaucoma. The implant consists of a medical grade silicone tube secured near the anterior suture points of the medical grade silicone episcleral plate. Two models are available that cover surface areas of approximately 250 mm² and 350 mm². The anterior suture points are located on the anterior side of the plate, flanking each side of the tubing track. The silicone plate is barium impregnated and contains a minimum of 8% and maximum of 10% of barium sulfate (BaSO₄) to increase ultrasound resolution and identification with CT scans, MRIs and plain skull films. The plate conforms to the shape of the globe at its equator and


provides a surface from which fluid can be dispersed. The device is for single use only.

Each Ahmed ClearPath® glaucoma drainage device is supplied with a 23-gauge hypodermic needle and a 4-inch 4-0 polypropylene *rip*cord, for optional use, in a sterile, sealed double-pouch. The use or non-use of both the 23-gauge hypodermic needle and/or the *rip*cord do not alter the standard care for the implantation of the Ahmed ClearPath® glaucoma drainage device. Once implanted, the Ahmed ClearPath® remains permanently in the patient's eye. The device is expected to be stable over the lifetime of the patient.

Model CP250

Model CP350

Model	Surface Area*	Plate Length	Plate Width (Max)	Plate Thickness (Max)	Tube Length	Tube ID	Tube OD
CP250	250 mm ²	16.70 mm	18.47 mm	0.86 mm	32 mm	0.305 mm	0.635 mm
CP350	350 mm ²	16.48 mm	30.47 mm	0.86 mm	32 mm	0.305 mm	0.635 mm

* Surface area does not include the area covered by the suture arms.

Indications for Use:

The Ahmed ClearPath® glaucoma drainage devices are indicated for the management of refractory glaucoma where previous surgical treatment has failed or is not expected to provide satisfactory

results. Such refractory glaucoma may include but is not limited to: neovascular glaucoma, congenital or infantile glaucoma, and refractory glaucoma resulting from aphakia or uveitis.

Intended Patient Population:

- Adult and pediatric patients

Clinical Benefits:

Intraocular pressure (IOP) reduction and/or reduction in the need for IOP-lowering medications. IOP is the only modifiable risk factor for glaucoma progression. Glaucoma progression can lead to blindness. The drainage tube is initially occluded by surgeons to restrict fluid outflow until scar tissue has formed around the implant. The scar tissue can provide sufficient resistance to fluid outflow, reducing the likelihood of hypotony. The benefits of the device are typically realized after the tube occlusion is removed. These benefits may diminish over time due to the natural wound-healing response around the implant which can restrict fluid outflow to the surrounding tissue. Consistent with the European Glaucoma Society (EGS) guidelines: "there is no single target IOP level that is appropriate for every patient, so the target IOP needs to be estimated separately for each eye of every patient and adapted to changing patient status". Further, the EGS states "Target IOP should be updated at each monitoring visit on the basis of changes in glaucoma or other ocular or systemic diseases." The document provides factors to consider when setting the target IOP. "There is little evidence to support any particular algorithm to set the target IOP. In newly diagnosed patients, the target IOP is initially determined according to the stage of disease and the baseline IOP. The treatment goal is typically estimated as a specific pressure level or a percentage reduction." It is up to the healthcare professional to judge the level of IOP sufficient to prevent the risk of progression of sight loss.

Summary of Safety and Clinical Performance:

The summary of safety and clinical performance (SSCP) for Ahmed ClearPath® is available in the European database on medical devices (EUDAMED) and it is linked to the basic UDI-DI: 00892064002348 (Model 250) and UDI-DI: 00892064002355 (Model 350). <https://ec.europa.eu/tools/eudamed>

Contraindications:

Bacterial conjunctivitis, bacterial corneal ulcers, endophthalmitis, no light perception, orbital cellulitis, bacteremia or septicemia and active scleritis.

Possible Complications and Adverse Events:

Serious or hemorrhagic choroidal detachments, hyphema, hypotony, flat or shallow anterior chamber, phthisis bulbi, retinal detachment, endophthalmitis, conjunctival erosion, wound leak, tube contacting cornea, tube block by iris or vitreous, tube migration anteriorly or posteriorly, corneal edema, need for secondary surgical intervention, uveitis and diplopia.

Warnings & Precautions:

1. Do not use the device if the sterile package integrity has been compromised.
2. Do not re-sterilize the implant.
3. Tampering with the implant can cause malfunctions.
4. Do not use the 4-inch 4-0 polypropylene *rip cord* as a suture.
5. The Ahmed ClearPath® glaucoma drainage devices are intended for single use only.
6. Use of surgical tools that can/may cause damage to the device are not recommended.
7. Store Ahmed ClearPath® at room temperature.

Storage Conditions:

Store Ahmed ClearPath® at room temperature.

Caution: Do not expose to direct sunlight and/or water

How Supplied/Expiration Date:

Each Ahmed ClearPath® glaucoma drainage device is supplied with a 23-gauge hypodermic needle and a 4-inch 4-0 polypropylene *rip cord* (preloaded in the tube lumen), for optional use, in a sterile, sealed double-pouch. Both pouches are held within a product specific cardboard box. Product Information, Patient Chart Labels, and Patient ID Card are also enclosed within the sterile package. The device has been terminally sterilized by gamma radiation. Sterility is assured provided that the peel pouch has not been opened or damaged and the sterility expiration date has not lapsed. The peel pouch is intended to be opened using a sterile technique, allowing the implant to be placed into the sterile field. The expiration date is indicated on the device's packaging. The implant should not be used after the expiration date indicated. **Note: The manufacturer disclaims all warranties expressed or implied, including but not limited to suitability for a particular purpose.**

Directions for Use:

Opening device and maintaining Sterile Field:

To open packaging and perform device transfer to the surgeon and sterile field, the following sequence of steps should be performed.

1. Verify the model number to be used and that the expiration date has not passed.
2. Open the carton, remove, and inspect the contents for damage using aseptic techniques.
3. Remove the inner pouch by peeling open the outer pouch using the flap adjacent to the chevron-shaped seal using aseptic techniques.

- Repeat the process for opening the inner pouch and remove the trays containing the devices using aseptic techniques.
 - Remove the tray cover and middle retainer securing the devices using aseptic techniques.
 - Present the tray to the surgeon for retrieval of the device and provided needle using aseptic techniques or invert the tray allowing them to fall onto the sterile field.
 - The surgeon should examine the implant prior to use for integrity and proper size.
 - Examine the device to ensure particles are not present, as it may become electrostatically charged upon opening of the package. Rinse the device in sterile saline if required.
- (Fig. 4) The drainage tube is trimmed to permit a 2-3 mm insertion of the tube into the anterior chamber. The tube should be bevel cut to an anterior angle of 30° to facilitate insertion.
 - (Fig. 5) The AC is entered at 2-3 mm away from the limbus with a sharp 23-gauge needle to create a needle track, parallel to the iris. In addition, a paracentesis should be performed to allow for quick AC reformation in case of a flattened anterior chamber.

Caution: Care must be taken to ensure that the drainage tube does not contact the iris or corneal endothelium after insertion.

- (Fig. 6) The drainage tube is inserted approximately 2-3 mm into the anterior chamber through the previously created needle track.
- After surgery, the body's wound healing mechanism results in the formation of a fibrous capsule (bleb) composed of scar tissue around the implant, this scar tissue resists fluid outflow. In approximately 4-6 weeks the 7-0 vicryl suture dissolves. At that time, under a magnified view of a slit lamp biomicroscope or operating microscope the conjunctiva can be entered over the temporal side of the *ripcord* and the *ripcord* be removed using jeweler's forceps. Following the *ripcord* removal, the IOP should be checked and slit lamp and fundoscopic examination should be performed.

Implantation:

The surgeon should have knowledge and experience in the implantation of glaucoma drainage devices, generally, and the Ahmed family of devices, specifically, prior to implantation of the Ahmed ClearPath® glaucoma drainage device.

- (Fig. 1 Optional *Ripcord*) The tube must be ligated near the tube-plate junction preferably with a 7-0 vicryl dissolvable suture to occlude the tube. One technique to occlude the tube involves the placement of a releasable *ripcord* in the tube lumen at the time of the initial surgery. A 7-0 vicryl dissolvable suture is placed around the tube (and consequently around the *ripcord*) and tightened to prevent any flow through the tube and around the *ripcord*. A 4-inch 4-0 polypropylene *ripcord* (preloaded in the tube lumen) is provided with the implant.
- (Fig. 2) A conjunctival peritomy incision is made. A pocket is formed superiorly, avoiding the superior rectus and oblique muscles, with blunt dissection of Tenon's capsule from the episclera.
- (Fig. 3) The implant is placed away from the limbus under the conjunctiva and tenon's capsule between the rectus muscles (the endplate of model CP350 is wider than the model CP250 endplate and has portions that fit under adjacent rectus muscles) and sutured to the episclera. The leading edge of the endplate should be ~ 8-10 mm from the limbus resulting in the suture arms being 6-8 mm from the limbus.

Explant:

If removal of the device is necessary, standard surgical technique can be used to isolate the plate and tube from surrounding conjunctiva and Tenon's capsule.

- Any adhesions between the device and surrounding tissues should be amputated.
- Pre-existing sutures are cut and removed to allow for free mobility of the plate and tube.
- The silicone tube is then grasped and removed from the sclerotomy tunnel.
- The entire device is then removed from the operative site.
- The sclerotomy tunnel is sutured closed and water tight followed by re-approximation of the conjunctiva and Tenon's capsule using the surgeon's preferred suture and technique.

Directions for Use: (continued)

- Address and website of the manufacturer
- Website for patient information

Disposal:

1. Used disposable needles are contaminated and can transmit diseases. Immediately after use, dispose of the needle in an appropriately labeled puncture resistant biohazard sharps container. Dispose of *ripcord* and/or explanted devices in properly labeled Biohazard waste containers.
2. The Ahmed ClearPath® is a single patient, single use device. Do not reuse or re-sterilize. Re-sterilized or reused instruments pose risks of infection and may result in damage to tissue.

Return of Damaged/ Unsterile Product:

Contact New World Medical (NWM) for product return policy. Return the device with proper identification and reason for return. Label the return as a biohazard.

Contact Number: US: +1 (909) 466-4304

Return Address: New World Medical
10763 Edison Court
Rancho Cucamonga
California 91730, USA

Patient Information:

It is recommended that each patient receive proper information regarding glaucoma implants and/or glaucoma drainage devices prior to the decision to implant is made. The implanting surgeon should describe to the patient all the risks involved with glaucoma drainage device surgery. At the end of surgery, the Implant Card should be provided to the patient. The surgeon should give the patient the implant card provided with the appropriate information filled in, and should advise the patient to keep the card in a safe place, e.g. his or her wallet, for future reference. The surgeon should advise the patient that this implant card contains important information related to the Ahmed ClearPath® implant and that the card should be shown to their current and future health care providers. The surgeon should inform patient the implant card contains the following information:

- Implant date
- Name and Address of the implanting healthcare institution / provider
- Serial number
- Lot number
- The Unique Device Identifier (UDI)
- Model number

The surgeon should advise the patient that the Ahmed ClearPath® is comprised of materials with a long history of safe clinical use that remain stable within the ocular environment over the lifetime of the patient. The surgeon should provide information to the patient that there is no device specific post-op care outside of normal generic ophthalmic instructions.

Reporting:

Adverse reactions and/or potentially sight-threatening complications that may be reasonably regarded as product related and that were not previously expected in nature, severity or incidence must be reported to NWM. Report any U.S. product related event to the toll free NWM number at +1 (800) 832-5327. For customers from countries outside the U.S. please contact your local company representative. For customers from countries within the EU, any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Disclaimer of Liability:

NWM is not responsible for any technique(s) performed which are not stated in the instructions for use to implant the product.

Clinical Data:

Although the Ahmed ClearPath® has not been clinically evaluated for performance, the following reference in the published clinical literature summarizes the performance of the predicate device BAERVELDT® Glaucoma Implant.

o Budenz, Donald L., et al. "Five-year treatment outcomes in the Ahmed Baerveldt comparison study." Ophthalmology 122.2 (2015): 308-316.

The following additional clinical literature is referenced in support of differences in design of the Ahmed ClearPath® as compared to the predicate BAERVELDT® Glaucoma Implant.

- Placement of Device: This study was conducted to determine optimal positioning of various glaucoma drainage devices posterior to the limbus while maintaining a safe distance relative to the optic nerve.

o Kahook, Malik Y., et al. "Location of glaucoma drainage devices relative to the optic nerve."

British Journal of Ophthalmology 90.8 (2006): 1010-1013.

- Implant Design – Ridge: The following studies evaluated the safety and effectiveness of the AGV-FP7 reference device to the predicate BAERVELDT® Glaucoma Implant and indicate that a lack of a ridge does not significantly impact implant performance.

o Budenz, Donald L., et al. "Five-year treatment outcomes in the Ahmed Baerveldt comparison study." *Ophthalmology* 122.2 (2015): 308-316.



o Budenz, Donald L., et al. "Postoperative complications in the Ahmed Baerveldt comparison study during five years of follow-up." *American Journal of Ophthalmology* 163 (2016): 75-82.

- Use of *Ripcord* for Tube Occlusion: A *ripcord* for optional use is provided within the tube lumen of the Ahmed ClearPath device. Tube occlusion via intraluminal *ripcord* is a common technique to prevent early postoperative hypotony in non-valved glaucoma drainage devices. Below are examples of medical publications that discuss this technique, its effectiveness and its impact on implant performance.

o Sherwood, Mark B., and M. Fran Smith. "Prevention of early hypotony associated with Molteno implants by a new occluding stent technique." *Ophthalmology* 100.1 (1993): 85-90.

o An, Selena J., et al. "Scheduled postoperative ripcord removal in Baerveldt 350 implants corresponds with increased complications: a prospective, randomized trial." *Investigative Ophthalmology & Visual Science* 59.9 (2018): 2073-2073.

Symbols Used on Packaging:

	Manufacturer		Authorized representative in the European Community		Sterilized using irradiation
	Date of manufacture		Consult instructions for use or consult electronic instructions for use		Do not use if package is damaged
	Use-by date		Read all warnings and precautions in instructions for use		MR safe
	Serial number		Do not reuse		Quantity of Device in Package
	Batch Code		Do not re-sterilize		Catalogue Number
	Single Sterile Barrier system with protective packaging outside		Distributor		Medical Device
	Double Sterile Barrier System		Keep Dry		Keep away from sunlight
	Single Sterile Barrier System with Protective packaging inside		Caution: Federal (USA) law restricts this device to sale by or on the order of a physician		Unique Device Identifier
	Open Here				

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o Budenz, Donald L., et al. Five-year treatment outcomes in the Ahmed Baerveldt comparison study. Ophthalmology 122.2 (2015): 308-316.

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o Kahook, Malik Y., et al. Location of glaucoma drainage devices relative to the optic nerve. British Journal of Ophthalmology 90.8 (2006): 1010-1013.

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o Budenz, Donald L., et al. Five-year treatment outcomes in the Ahmed Baerveldt comparison study. *Ophthalmology* 122.2 (2015): 308-316.



o Budenz, Donald L., et al. Postoperative complications in the Ahmed Baerveldt comparison study during five years of follow-up. *American Journal of Ophthalmology* 163 (2016): 75-82.

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o Sherwood, Mark B., and M. Fran Smith. Prevention of early hypotony associated with Molteno implants by a new occluding stent technique. *Ophthalmology* 100.1 (1993): 85-90.

o An, Selena J., et al. Scheduled postoperative ripcord removal in Baerveldt 350 implants corresponds with increased complications: a prospective, randomized trial. *Investigative Ophthalmology & Visual Science* 59.9 (2018): 2073-2073.

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California, 91730 U.S.A.
Emergo Europe B.V.