# The Ahmed<sup>®</sup> Glaucoma Valve

Model FP8

NEW WORLD  
MEDICAL

CE  
0459

# Index

|                            |     |
|----------------------------|-----|
| English (en) . . . . .     | 3   |
| Français (fr) . . . . .    | 9   |
| Deutsch (de) . . . . .     | 15  |
| Italiano (it) . . . . .    | 21  |
| Español (es) . . . . .     | 27  |
| Português (pt) . . . . .   | 33  |
| Čeština (cs) . . . . .     | 39  |
| Dansk (da) . . . . .       | 45  |
| Ελληνικά (el) . . . . .    | 51  |
| Suomi (fi) . . . . .       | 57  |
| Hrvatski (hr) . . . . .    | 63  |
| Latviski (lv) . . . . .    | 69  |
| Magyar (hu) . . . . .      | 75  |
| Nederlands (nl) . . . . .  | 81  |
| Polski (pl) . . . . .      | 87  |
| Română (ro) . . . . .      | 93  |
| Slovenščina (sl) . . . . . | 99  |
| Svenska (sv) . . . . .     | 105 |
| Lietuvių (lt) . . . . .    | 111 |
| Norsk (no) . . . . .       | 117 |
| Íslenska (is) . . . . .    | 123 |

English (en)

# The Ahmed® Glaucoma Valve Model FP8

## Instructions for Use

### DESCRIPTION:

The Ahmed® Glaucoma Valve Model FP8 (AGV-FP8) is a valved aqueous drainage implant designed to regulate intraocular pressure in eyes suffering from refractory glaucoma. The Ahmed® device consists of a valve mechanism and an end plate. The valve mechanism is comprised of a silicone drainage tube and polypropylene casing which houses a silicone elastomer valve membrane. The polypropylene casing protects the valve membrane from blockage by fibrous tissue. The end plate, made from silicone, conforms to the shape of the globe at its equator and provides a surface from which fluid can be dispersed. Due to a smaller end plate surface area compared to other glaucoma drainage devices such as the AGV-FP7, the AGV-FP8 can be used in situations where a smaller end plate device is desired. The device is for single use only. Once implanted, the Ahmed® Glaucoma Valve remains permanently in the patient's eye. The device is expected to be stable over the lifetime of the patient.

### INDICATIONS:

The Ahmed® Glaucoma Valve Model FP8 is indicated for the management of refractory glaucomas, where previous surgical treatment has failed, or by experience is known not to provide satisfactory results. Such refractory glaucomas can include neovascular glaucoma, primary open angle glaucoma unresponsive to medication, congenital or infantile glaucoma, and refractory glaucomas resulting from aphakia or uveitis.

### INTENDED PATIENT POPULATION:

- Adult and pediatric patients.

### CLINICAL BENEFITS:

Immediate reduction in intraocular pressure (IOP) and/or reduction in the need for IOP lowering medications. IOP is the only modifiable risk factor for glaucoma progression. Glaucoma progression can lead to blindness. These clinical benefits may diminish over time due to the natural wound healing response around the implant which restricts fluid outflow to the surrounding tissue. The valve mechanism reduces the likelihood of hypotony compared to non-valved glaucoma drainage devices which do not have a built-in outflow restriction mechanism.

Consistent with the European Glaucoma Society (EGS) guidelines: "there is no single target IOP level that is appropriate

for every patient, so the target IOP needs to be estimated separately for each eye of every patient and adapted to changing patient status". Further, the EGS states "Target IOP should be updated at each monitoring visit on the basis of changes in glaucoma or other ocular or systemic diseases."

The document provides factors to consider when setting the target IOP. "There is little evidence to support any particular algorithm to set the target IOP. In newly diagnosed patients, the target IOP is initially determined according to the stage of disease and the baseline IOP. The treatment goal is typically estimated as a specific pressure level or a percentage reduction." It is up to the healthcare professional to judge the level of IOP sufficient to prevent the risk of progression of sight loss.

### SUMMARY OF SAFETY AND CLINICAL PERFORMANCE:

The summary of safety and clinical performance (SSCP) for Ahmed Glaucoma Valve is available in the European database of medical devices (EUDAMED) and it is linked to the basic UDI-DI: 00892064002126.

<https://ec.europa.eu/tools/eudamed>

### CONTRAINDICATIONS:

- |                             |                             |
|-----------------------------|-----------------------------|
| a. Bacterial conjunctivitis | e. Bacteremia or septicemia |
| b. Bacterial corneal ulcers | f. Active scleritis         |
| c. Endophthalmitis          | g. No light perception      |
| d. Orbital Cellulitis       |                             |

### COMPLICATIONS AND ADVERSE REACTIONS:

Complications and adverse reactions during or following surgery may include:

- |                      |                      |
|----------------------|----------------------|
| Corneal Edema        | Choroidal detachment |
| Corneal Touch        | Iritis               |
| Iris/Tube Touch      | Hyphema              |
| Synechia             | Tube Obstruction     |
| Exposed Sclera Graft | Tube Retraction      |

as well as known complications of aqueous shunts and general intraocular surgery including hypotony,

shallow chamber, choroidal effusion, suprachoroidal hemorrhage, retinal detachment, cataract, conjunctival buttonhole, phthisis bulbi, bullous keratopathy, uveitis and endophthalmitis.

### **WARNINGS, PRECAUTIONS:**

Do not use the device if sterile package integrity has been compromised. Do not re-sterilize the implant. Before using the Ahmed® Glaucoma Valve, the implanting surgeon should be skilled in glaucoma filtering procedures and familiar with the use of drainage devices, as well as post-operative care required.

Priming the valve is essential for proper functioning of the valve. Do not implant the valve without priming it.

Ensure that the valve is primed before operating on the patient. Inherently, the force required to prime the valve may vary between units. Some valves may require more force applied over a longer period of time than other valves in order to perform the priming operation. If you are unable to prime the valve, do not use the device. The manufacturer will exchange the device.

In rare cases, you may be unable to prime the valve due to crimped or pinched tubing. In such a situation, do not use the device. The manufacturer will exchange the device.

Using a needle for priming can puncture the tube resulting in an undesirable leak or inability to prime. Only use a blunt cannula to prime the valve.

The Ahmed® Glaucoma Valve is intended for single use, it shall not be reused. It shall be discarded to avoid deterioration, cross contamination and/or infection.

Tampering with the valve can cause malfunctioning of the valve.

### **PATIENT INFORMATION:**

The surgeon should give the patient the implant card provided with the appropriate information filled in, and should advise the patient to keep the card in a safe place, e.g., his or her wallet, for future reference. The surgeon should advise the patient that this implant card contains important information related to the Ahmed® Glaucoma Valve implant and that the card should be shown to their current and future health care providers. The surgeon should inform patient the implant card contains the following information:

- Implant date
- Name and Address of the implanting healthcare institution / provider
- Serial number
- Lot number

- The Unique Device Identifier (UDI)
- Model number
- Address and website of the manufacturer
- Website for patient information

The surgeon should advise the patient that the Ahmed® Glaucoma Valve is comprised of materials with a long history of safe clinical use that remain stable within the ocular environment over the lifetime of the patient. The surgeon should provide information to the patient that there is no device specific post-op care outside of normal generic ophthalmic instructions.

### **REPORTING:**

Adverse reactions and/or potentially sight-threatening complications that may be reasonably regarded as product related and that were not previously expected in nature, severity or incidence must be reported to New World Medical (NWM). Report any U.S. product related event to the toll free NWM number at +1 (800) 832-5327. For customers from countries outside of the U.S. please contact your local company representative. For countries within the EU, any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

### **STORAGE CONDITIONS:**


Store Ahmed® Glaucoma Valve at room temperature.

**Caution:** Do not expose to direct sunlight and/or water.


## AGV-FP8

Trade Name: Ahmed® Glaucoma Valve Model FP8

### Top View

### Side View

**End Plate:** Medical Grade Silicone

**Drainage Tube:** Medical Grade Silicone

**Valve Membrane:** Medical Grade Silicone

**Valve Casing:** Medical Grade Polypropylene (PP)

**Adhesive:** Medical Grade Silicone Glue

**Maximum Thickness:** 2.1 mm

**Max Length:** 36.4 mm

**Plate Length:** 11.0 mm

**Device Overall Width:** 11.0 mm

**Drainage Area/Explant Surface Area:** 102 mm<sup>2</sup>

The implant is supplied sterile in a sealed pouch. Product Information, Patient Chart Labels, Implant Card, and Implant Card Instructions are also enclosed with the sterile package. The implant has been terminally sterilized by Gamma Radiation. Sterility is assured provided that the peel pouch has not been opened or damaged, and the sterility expiration date has not lapsed. The peel pouch is intended to be opened using sterile technique, allowing the implant to be dropped into the sterile field.

#### OPENING DEVICE AND MAINTAINING STERILE FIELD:

To open packaging and perform device transfer to the surgeon and sterile field, the following sequence of steps should be performed.

- Verify the model number to be used and that the expiration date has not passed.
- Open the carton, remove, and inspect the contents for damage using aseptic techniques.
- Remove the inner pouch by peeling open the outer pouch using the flap adjacent to the chevron-shaped seal using aseptic techniques.


- Repeat the process for opening the inner pouch and remove the implant using aseptic techniques.
- The surgeon should examine the implant prior to use for integrity and to ensure particles are not present, as it may become electrostatically charged upon opening of the package. Rinse the device in sterile saline if required.

#### DISPOSAL:

1. Dispose of explanted devices in properly labeled Biohazard waste containers.
2. The Ahmed® Glaucoma Valve is a single patient, single use device. Do not reuse or re-sterilize. Resterilized or reused instruments pose risks of infection and may result in damage to tissue.

## SURGICAL PROCEDURE:

### Figure 1

Refer to figure 1 for the illustration of each step. The steps described here are intended as a guideline only, and do not represent recommended treatment for any particular patient. The use of a specific surgical technique or maneuver is at the sole discretion of the surgeon. Surgeons should be familiar with the use of glaucoma drainage devices and post-operative care considerations before implanting a drainage device. Reference papers and surgical video tapes are available upon request.

1. The implant should be examined and primed prior to implantation. Priming is accomplished by injecting 1 mL of balanced salt solution or sterile water through the drainage tube and valve, using a blunt 26-gauge cannula.
2. A conjunctival peritomy incision is made. A pocket is formed superiorly, avoiding the superior rectus and oblique muscles, with blunt dissection of Tenon's capsule from the episclera.
3. The valve body is inserted into the pocket between the rectus muscles and sutured to the episclera. Rotate the suture knots into the the suture hole openings. The leading edge of the end plate should be ~ 8-10 mm from the limbus.
4. The drainage tube is trimmed to permit a 2-3 mm insertion of the tube into the anterior chamber (AC).

The tube should be bevel cut to an anterior angle of 30° to facilitate insertion.

5. The AC is entered at 1-2 mm way from the limbus with a sharp 23-gauge needle to create a needle track, parallel to the iris. In addition, a paracentesis should be performed to allow for quick AC reformation in case of a flattened AC. Caution: Care must be taken to ensure that the drainage tube does not contact the iris or corneal endothelium after insertion.

Note: Some surgeons prefer to enter the AC from at least 3 mm away from the limbus.

6. The drainage tube is inserted approximately 2-3 mm into the AC through the previously created needle track. The leading edge of the end plate should be 8-10 mm from the limbus.
7. The exposed drainage tube is covered with a piece of preserved, donor sclera, pericardium, cornea, or other suitable patch graft material which is sutured into place and the conjunctiva is closed.

Note: As an alternative to the use of patch graft material, a 2/3 thickness limbal-based scleral flap may be made. The tube is inserted into the AC through a 23-gauge needle puncture made under the flap. The flap is then sutured closed.

## EXPLANT:

If removal of the device is necessary, standard surgical technique can be used to isolate the plate and tube from surrounding conjunctiva and Tenon's capsule.

1. Any adhesions between the device and surrounding tissues should be amputated.
2. Pre-existing sutures are cut and removed to allow for free mobility of the plate and tube.
3. The silicone tube is then grasped and removed from the sclerotomy tunnel.
4. The entire device is then removed from the operative site.
5. The sclerotomy tunnel is sutured closed and watertight followed by re-approximation of the conjunctiva and Tenon's capsule using the surgeon's preferred suture and technique.

## CLINICAL STUDY RESULTS:

Published clinical data for the AGV-FP8 can be found in the citations referenced below. This includes safety and effectiveness data and information on adverse effects and complications associated with the AGV-FP8.

1. Comparison of the Outcome of Silicone Ahmed Glaucoma Valve Implantation with a Surface Area between 96 and 184 mm<sup>2</sup> in Adult Eyes. This study was a retrospective review of records from adult refractory glaucoma patients who underwent either AGV-FP8 or AGV-FP7 implantation by two surgeons at a single center. Similar surgical techniques were used regardless of implant type. Some patients were followed up to 3 years after surgery. There were no statistically significant differences between the groups in preservation of vision, IOP reduction, or decrease in the number of glaucoma medications.

Koh KM, Hwang YH, Jung JJ, Sohn YH, Kim HK. Comparison of the outcome of silicone Ahmed glaucoma valve implantation with a surface area between 96 and 184 mm<sup>2</sup> in adult eyes. *Korean J Ophthalmol*. 2013 Oct;27(5):361-7. PubMed PMID: 24082774; PubMed Central PMCID: PMC3782582.

2. Outcomes of Ahmed Valve Implant Following a Failed Initial Trabeculectomy and Trabeculectomy in Refractory Primary Congenital Glaucoma. This was a retrospective noncomparative case series of eyes with a diagnosis of refractory primary congenital glaucoma. The AGV-FP8 was implanted by a single surgeon after a failed primary trabeculectomy + trabeculectomy. In this difficult-to-treat group of patients, the AGV-FP8 was found to be an effective treatment.

Dave P, Senthil S, Choudhari N, Sekhar GC. Outcomes of Ahmed valve implant following a failed initial trabeculectomy and trabeculectomy in refractory primary congenital glaucoma. *Middle East Afr J Ophthalmol*. 2015 Jan-Mar;22(1):64-8. PubMed PMID: 25624676; PubMed Central PMCID: PMC4302479.

3. Combined trabeculectomy-trabeculectomy versus Ahmed valve implantation for refractory primary congenital glaucoma in Egyptian patients: a long-term follow-up. This was a randomized, prospective, single-surgeon, comparative study that included 66 eyes with refractory primary congenital glaucoma with up to four years of follow up reported. Patients had previously failed goniotomy and trabeculectomy. Half of the patients underwent a combined trabeculectomy-trabeculectomy procedure and the other half underwent AGV-FP8 implantation. Both procedures were found to be suitable options in advanced refractory primary congenital glaucoma with similar long term IOP reduction, decrease in number of glaucoma medications, and success rates. A higher rate of hyphema was reported in the combined trabeculectomy-trabeculectomy group with other rates of complications similar between the two groups.

Helmy, Hazem. "Combined trabeculectomy-trabeculectomy versus Ahmed valve implantation for refractory primary congenital glaucoma in Egyptian patients: a long-term follow-up." *Electronic physician* 8.2 (2016): 1884.

4. Surgical outcomes of additional Ahmed glaucoma valve implantation in refractory glaucoma. Clinical histories of 23 refractory glaucoma patients, 21 of whom underwent a AGV-FP8 implantation after a failed glaucoma drainage device implantation were retrospectively reviewed. Outcomes for up to 3 years were reported. Implantation of an AGV-FP8 was described as a good choice for surgical treatment when the first glaucoma drainage device failed. Corneal decompensation was found in some cases but no other serious complications were reported.

Ko, Sung Ju, et al. "Surgical outcomes of additional Ahmed glaucoma valve implantation in refractory glaucoma." *Journal of glaucoma* 25.6 (2016): e620-e624.

## RETURN OF DAMAGED/UNSTERILE PRODUCT:

Contact NWM for product return policy. Return the device with proper identification and reason for return. Label the return as a biohazard.

Contact Number: US: +1 (909) 466-4304

Email: [orders@newworldmedical.com](mailto:orders@newworldmedical.com)

Return Address: New World Medical  
10763 Edison Court  
Rancho Cucamonga  
California 91730, USA

## Symbols Used on Packaging:

|  |   |   |  |   |                                  |
|--|---|---|--|---|----------------------------------|
|  | Manufacturer  |  | Authorized representative in the European Community  |  | Sterilized using irradiation     |
|  | Date of manufacture   |  | Consult instructions for use or consult electronic instructions for use                    |  | Do not use if package is damaged |
|  | Use-by date   |  | Read all warnings and precautions in instructions for use                                  |  | MR safe                          |
|  | Serial number   |  | Do not reuse   |  | Quantity of Device in Package    |
|  | Batch Code  |  | Do not re-sterilize  |  | Catalogue Number                 |
|  | Single Sterile Barrier system with protective packaging outside |  | Distributor  |  | Medical Device                   |
|  | Double Sterile Barrier System                                   |  | Keep Dry   |  | Keep away from sunlight          |
|  | Single Sterile Barrier System with Protective packaging inside  |  | Caution: Federal (USA) law restricts this device to sale by or on the order of a physician |  | Unique Device Identifier         |
|  | Open Here   |   |  |   |                                  |

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Koh KM, Hwang YH, Jung JJ, Sohn YH, Kim HK. Comparison of the outcome of silicone Ahmed glaucoma valve implantation with a surface area between 96 and 184 mm<sup>2</sup> in adult eyes. Korean J Ophthalmol. 2013 Oct;27(5):361-7. PubMed PMID: 24082774; PubMed Central PMCID: PMC3782582.

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Helmy, Hazem. „Combined trabeculectomy-trabeculectomy versus Ahmed valve implantation for refractory primary congenital glaucoma in Egyptian patients: a long-term follow-up.“ *Electronic physician* 8.2 (2016): 1884.

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Koh KM, Hwang YH, Jung JJ, Sohn YH, Kim HK. Comparison of the outcome of silicone Ahmed glaucoma valve implantation with a surface area between 96 and 184 mm<sup>2</sup> in adult eyes. *Korean J Ophthalmol.* 2013 Oct; 27(5): 361-7. PubMed PMID: 24082774; PubMed Central PMCID: PMC3782582.

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Helmy, Hazem. "Combined trabeculotomy-trabeculectomy versus Ahmed valve implantation for refractory primary congenital glaucoma in Egyptian patients: a long-term follow-up." *Electronic physician* 8.2 (2016): 1884.

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Koh KM, Hwang YH, Jung JJ, Sohn YH, Kim HK. Comparison of the outcome of silicone Ahmed glaucoma valve implantation with a surface area between 96 and 184 mm<sup>2</sup> in adult eyes. *Korean J Ophthalmol.* 2013 okt;27(5):361-367. PubMed PMID: 24082774; PubMed Central PMCID: PMC3782582.

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Helmy, Hazem. 'Combined trabeculotomy-trabeculectomy versus Ahmed valve implantation for refractory primary congenital glaucoma in Egyptian patients: a long-term follow-up.' *Electronic physician* 8.2 (2016): 1884.

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Ko, Sung Ju, et al. 'Surgical outcomes of additional Ahmed glaucoma valve implantation in refractory glaucoma.' *Journal of glaucoma* 25.6 (2016): e620-e624.

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California, 91730 U.S.A.
Emergo Europe B.V.