KDB GLIDE®

## Instructions For Use

CE  
0459

# **Index**

English (en) . . . . .	1
Čeština (cs) . . . . .	11
Dansk (da) . . . . .	22
Deutsch (de) . . . . .	33
Ελληνικά (el) . . . . .	44
Español (es) . . . . .	56
Suomi (fi) . . . . .	67
Français (fr) . . . . .	78
Hrvatski (hr) . . . . .	89
Magyar (hu) . . . . .	100
Íslenska (is) . . . . .	111
Italiano (it) . . . . .	122
Lietuvių k. (lt) . . . . .	133
Latviski (lv) . . . . .	143
Nederlands (nl) . . . . .	153
Norsk (no) . . . . .	164
Polski (pl) . . . . .	175
Português (pt) . . . . .	186
Română (ro) . . . . .	197
Slovenščina (sl) . . . . .	208
Svenska (sv) . . . . .	219

	English (en)
<b>A</b>	Shaft
<b>B</b>	Indicator for Blade Orientation
<b>C</b>	Dual Blades
<b>D</b>	Pointed Tip
<b>E</b>	Heel
<b>F</b>	Ramp

 **CAUTION**

- A. Physician training is required before use of this instrument.
- B. Contents are guaranteed sterile provided that the package including the tray seal is not punctured or damaged. Do not use if the package is damaged or open, or if the seal appears damaged.
- C. Read instructions before use.
- D. Do not use if the product is past its expiration date.
- E. When removing the instrument from its tray ensure that the tip does not touch any part of the tray. The blade tip may dull if it touches the tray surface. Do not touch the tip to any surfaces other than those of the eye as it may damage or dull the tip.
- F. Do not use the instrument if it contains debris, is dull, bent, rusted or otherwise damaged.
- G. Precautions standard to ophthalmic surgery shall be observed when using the KDB GLIDE®, including customary aseptic technique.
- H. The instrument contains stainless steel. Do not use for patients sensitive or allergic to metals.
- I. Used sharps are contaminated and can transmit diseases. Immediately after use, dispose of the instrument in an appropriately labeled puncture resistant biohazard sharps container.
- J. The KDB GLIDE® instrument is a single patient, single use disposable instrument. Do not reuse or re-sterilize. Re-sterilized or reused instruments pose risks of infection and dull blades that may result in damage to tissue.

**DESCRIPTION AND INTENDED USE**

The KDB GLIDE® is a disposable ophthalmic knife, a manual surgical instrument used in ophthalmic surgical procedures to surgically remove a strip of trabecular meshwork (TM) in adult patients. The instrument is sterile and is intended for single patient, single use.

The KDB GLIDE® consists of a surgical grade stainless steel body and a plastic handle. The stainless steel body is comprised of a long shaft that allows for access to the TM, a pointed tip to pierce TM, a ramp that lifts and stretches tissue as the device is advanced, and two blades that surgically remove a strip of TM. An indicator on the plastic handle specifies the direction to advance the instrument.

**INDICATIONS FOR USE**

The KDB GLIDE is indicated to surgically remove a strip of trabecular meshwork to reduce intraocular pressure in adult patients with primary open angle glaucoma.

**INTENDED POPULATION**

Adult patients.

**CONTRAINDICATIONS**

Like most angle surgeries, proper visualization using a gonioprism is key to surgically remove a strip of TM with the instrument. Do not attempt to use the instrument if corneal clarity is poor or if visualization of angle structures is not possible.

**NOTICE:** Any serious injury that has occurred in relation to this device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

**CAUTION:** Do not use KDB GLIDE® if you are not trained for intraoperative gonioscopy or have poor visualization of angle structures. Improper visualization may result in misuse of the instrument and damage to eye structures adjacent to the TM.

**UNDESIRABLE EFFECTS / SIDE EFFECTS**

Adverse events that may be reasonably associated with the use of KDB GLIDE® in the eye include but are not limited to the following: elevation of intraocular pressure, hyphema, Descemet's membrane tear or detachment, iridodialysis, iris

tear, iris injury, iritis, change to pupil shape, cyclodialysis cleft, corneal decompensation, corneal injury, corneal edema or opacification, perforation of sclera, as well as known complications of general intraocular surgery including anterior chamber shallowing, hypotony, hypotony maculopathy, IOL dislocation, cataract formation, posterior capsule opacification, loss of vitreous, choroidal detachment, choroidal effusion, suprachoroidal hemorrhage, vitreous hemorrhage, posterior vitreous detachment, retinal detachment, retinal dialysis, retinal tear, endophthalmitis or other ocular infection.

The following warnings and precautions are measures to be taken in the event of malfunction of the device in its performance that may affect safety.

## **WARNINGS**

- Do not use KDB GLIDE® if you are not trained for intraoperative gonioscopy
- The KDB GLIDE® is intended to surgically remove a strip of TM and cannot be used to make a clear corneal incision.

## **PRECAUTIONS**

- Contents are provided sterile and intended for single patient, single use. Do not re-sterilize or re-use the device.
- Re-sterilized or reused instruments pose risks of infection and dull blades that may result in damage to tissue.
- Do not use if package, including the tray seal, are punctured or damaged. The KDB GLIDE® is intended to surgically remove a strip of TM and cannot be used to make a clear corneal incision.
- Do not use past expiration date.
- When removing the instrument from its tray ensure that the tip does not touch any part of the tray. The blade tip may dull if it touches the tray surface. Do not touch the tip to any surfaces other than those of the eye as it may damage or dull the tip. Do not use KDB GLIDE® if you are not trained for intraoperative gonioscopy
- Do not use instrument if it contains debris, is dull, bent, rusted, or otherwise damaged.
- Observe all usual precautions undertaken during intraocular surgery.
- Use aseptic technique and ensure KDB GLIDE® and field sterility as is customary during intraocular surgery.
- The instrument contains stainless steel. Do not use for patients sensitive or allergic to metals.
- Used sharps are contaminated and can transmit diseases. Immediately after use, dispose of the instrument in an appropriately labeled puncture-resistant biohazard sharps container.

## **INSTRUCTIONS**

### **1. Removing the Instrument from the Plastic Tray**

- 1.1. **CAUTION:** See G.
- 1.2. Open the box and peel away the sealed Tyvek lid from the plastic container.
- 1.3. Remove the KDB GLIDE® from the plastic container taking care to lift the tip away from the tray first by gripping the KDB GLIDE® handle at the finger slots provided in the tray.
- 1.4. **CAUTION:** See E.
- 1.5. Inspect instrument. **CAUTION:** See F.

### **2. Considerations**

- 2.1. Like other angle surgeries, use of a gonioprism over the patient's cornea under surgical microscope magnification is required for visualization of the TM and the instrument during use.
- 2.2. Positioning of the patient is essential to obtaining a proper view.
- 2.3. The KDB GLIDE® can be used stand-alone or combined with other procedures.
- 2.4. The KDB GLIDE® is inserted into the eye through a previously created clear corneal incision. Larger incisions allow for greater maneuverability.  
**CAUTION:** The KDB GLIDE® is intended to surgically remove a strip of TM and cannot be used to make a clear corneal incision.
- 2.5. The clear corneal incision should be approximately 180 degrees away from the TM intended to be surgically removed.

## 4.

2.6. The anterior chamber should be inflated with viscoelastic when using the KDB GLIDE®. **CAUTION:** Overinflating the chamber may make cutting of the TM difficult as the Canal of Schlemm may collapse.

### 3. Using KDB GLIDE® under Gonioscopic View to surgically remove a strip of TM.

3.1. Allowing some blood to reflux into the Canal of Schlemm, by lowering intraocular pressure, can be used to identify TM. Insert the tip and shaft of KDB GLIDE® through the clear corneal incision and advance to the opposite side of the anterior chamber. **CAUTION:** Be careful to avoid touching the corneal endothelium and iris with the instrument as you are advancing or retracting the instrument in the anterior chamber. When inserting the tip, ensure that the long edge of the footplate is parallel to the incision (sharp tip is not perpendicular to incision slit). (**Figure a**)

3.2 Engage the TM across from the clear corneal incision with the pointed tip of the KDB GLIDE® and pierce through TM and into the Canal of Schlemm while avoiding the anterior wall. (**Figure b**)

3.3 Advance KDB GLIDE® three to five clock hours in the direction of the TM region that is intended to be surgically removed by pivoting and advancing/retracting the shaft along the clear corneal incision while keeping the KDB GLIDE® footplate in the Canal of Schlemm. (**Figure c**) Avoid cutting surrounding tissue. Some blood reflux into the anterior chamber from the collector channels can be expected as the TM is removed. Viscoelastic or BSS can be used to improve visualization of the area being treated. Blood can be aspirated out of the anterior chamber. **CAUTION:** Ensure that the smooth heel of the tip slides forward with minimal pressure against the anterior wall of the Canal of Schlemm. If the heel is angled away from the back wall, the pointed tip may engage the anterior wall of the Canal of Schlemm and prevent smooth advancement of the KDB GLIDE® or result in sub-optimal engagement of the targeted TM tissue.

3.4 Detach the TM strip using KDB GLIDE® or microforceps. To detach the TM strip with the KDB GLIDE®, engage the TM across from the clear corneal incision with the pointed tip and pierce through the TM and into the Canal of Schlemm while avoiding the anterior wall of the Canal of Schlemm. Advance KDB GLIDE® in one direction until the desired endpoint of goniotomy is reached. Reverse the blade and advance towards the original area of treatment until a connection is made with the initial area of goniotomy.

3.5 Do not leave the detached TM strip in the anterior chamber. The detached TM tissue can be removed using KDB GLIDE® as the tissue may adhere to the blade, using ophthalmic forceps, or through aspiration from the clear corneal wound. Only aspirate fully detached TM strips.

3.6 **CAUTION:** If the angle of approach of the KDB GLIDE® instrument to TM is too steep only one of the dual blades may incise TM.

3.7 After surgically removing a strip of TM, retract the KDB GLIDE® through the original clear corneal incision.

### 4. Disposal

4.1. **CAUTION:** See I.

4.2. **CAUTION:** See J.

### HOW SUPPLIED

One KDB GLIDE® is supplied per box and comes packaged in a plastic tray with a sealed Tyvek lid. The tray secures and protects the delicate tip of the KDB GLIDE®. The KDB GLIDE® is a disposable, single patient, single use, sterile surgical instrument.

KDB GLIDE® is sterilized with gamma radiation.

### STORAGE REQUIREMENTS

KDB GLIDE® should be stored at room temperature. **CAUTION:** Do not expose to water, sunlight, extreme temperatures, or high humidity.

### EXPIRATION DATE

The KDB GLIDE® can be used up to 3 years after the date of manufacture as indicated by the expiration date on the label. Sterility is guaranteed until the expiration date if the packaging is not punctured or otherwise damaged. **CAUTION:** Do not use the KDB GLIDE® past the indicated expiration date.

## CLINICAL DATA SUMMARY

The reference device, OMNI Surgical System (K202678), was used to support the scientific methodology to evaluate the ability of KDB to reduce IOP in adult patients with POAG and a systematic literature review was performed to characterize the safety of the predicate device, NMX-1000

### KDB Glide® Procedure Performed with Cataract Surgery (Falkenberry et al)

Falkenberry S, Singh IP, Crane CJ, Haider MA, Morgan MG, Grenier CP et al. Excisional goniotomy vs trabecular microbypass stent implantation: a prospective randomized clinical trial in eyes with mild to moderate open-angle glaucoma. J Cataract Refract Surg. 2020; 46(8):1165-71

Falkenberry et al conducted a prospective, randomized, active-controlled, parallel group clinical trial of the KDB Glide® procedure combined with cataract surgery. The study was conducted at 9 study sites in the U.S. and surgeries were performed between June 2016 and January 2019. The purpose of the study was to compare reduction in intraocular pressure (IOP) and IOP-lowering medication in eyes undergoing excisional goniotomy with Kahook Dual Blade (KDB) Glide® vs iStent trabecular micro-bypass (Glaukos, Inc.) implantation, both combined with phacoemulsification, in eyes with mild-to-moderate open-angle glaucoma (OAG).

The study consisted of eyes diagnosed with open-angle glaucoma, including pseudoexfoliative glaucoma and pigmentary glaucoma. The KDB-phaco cohort consisted of 82 eyes of 82 subjects. The majority of study eyes (n=75, 95.1%) had a diagnosis of primary open-angle glaucoma (POAG); 5 eyes (6.1%) had pseudoexfoliative glaucoma (PEG) and the remainder (n=2, 2.4%) had pigmentary glaucoma (PG). The glaucoma severity was characterized as mild (79.3%) to moderate (20.7%). Subjects were followed for 12 months and were evaluated 1 day, 1 week, and 1, 3, 6, and 12 months postoperatively.

Baseline characteristics were as follows:

Baseline characteristics	KDB-Phaco (n=82)
Mean Age ± SD	70.2 ± 0.91 years
Sex	Female: 44 (53.7%) Male: 38 (46.3%)
Ethnicity	White: 66 (80.5%) African American: 10 (12.2%) Asian: 3 (3.7%) Hispanic: 1 (1.2%) Other: 2 (2.4%)
Medicated IOP Mean ± SD	18.5 ± 0.36 mmHg
Glaucoma Medications Mean ± SD	1.31 ± 0.07 meds

The key inclusion criteria consisted of subjects 18-90 years of age with mild-to-moderate POAG, PEG, or PG who were on 1-3 topical ocular hypotensive medications. In order to qualify, subjects needed a baseline IOP of 14-28 mmHg (inclusive) and a visually significant cataract planned for elective extraction. Subjects using oral medications that could affect IOP, previous glaucoma surgery, recent ( $\leq$  3 months) glaucoma laser therapy, closed angles, or a history of steroid response were identified as key exclusion criteria. Under ophthalmic viscosurgical device and using a goniotrism for visualization, the KDB Glide® was introduced through a clear corneal incision approximately 180° from the targeted excision site. The pointed tip of the device engaged trabecular meshwork (TM), with the heel positioned within Schlemm's canal and the blade was advanced to excise a strip of trabecular meshwork (TM) approximately 3-4 clock hours in length using the parallel blades of the device. The KDB Glide® was then removed through the corneal incision.

The primary effectiveness endpoint was the proportion of eyes at 12 months with IOP reduction of 20% or greater or IOP medication reduction of 1 or more compared with baseline. Secondary effectiveness endpoints included the percentage reduction in IOP and the mean decrease in the number of IOP-lowering medications.

## 6.

Of the 82 eyes treated with KDB-phaco, 79 eyes were available for the 12-month visit and 3 eyes were considered lost to follow-up. There were no eyes with glaucoma secondary surgical intervention (SSI) for IOP control during the 12-month study; therefore, the effectiveness analyses included all 79 available eyes. The proportion of eyes at 12 months with an IOP reduction of 20% or greater from baseline or IOP medication reduction of 1 or more medication from baseline was 93.7% (74/79 eyes). Mean (SD) IOP decreased from a baseline of 18.5 ( $\pm$  0.4) mmHg to 15.4 ( $\pm$  0.4) mmHg at 12 months. Mean (SD) IOP-lowering medications were reduced from 1.3 ( $\pm$  0.1) at baseline to 0.3 ( $\pm$  0.1) at 12 months.

Safety outcomes included corrected distance visual acuity loss by 2 lines or more from best recorded to final visit, operative complications, device malfunctions, IOP elevations of 10 mmHg or greater from baseline, and adverse events. The following adverse events were reported:

Adverse Event	KDB-Phaco (n=82)
Increased IOP	26 eyes (31.7%)
Posterior Capsule Opacification	7 eyes (8.5%)
Hyphema	3 eyes (3.7%)
Cyclodialysis	1 eye (1.2%)

The most common adverse event was increased IOP, reported as occurring at least one postoperative visit in 26/82 (31.7%) eyes; these all resolved spontaneously or with medical management. Posterior capsule opacification (PCO) occurred in 7/82 (8.5%) eyes. Blood persisting in the anterior chamber beyond the first postoperative week (hyphema) occurred in 3/82 (3.7%) of study eyes. One eye (1.2%) developed a cyclodialysis cleft.

### KDB Glide® Procedure Performed Standalone (ElMallah et al)

ElMallah MK, Berdahl JP, Williamson BK, Dorairaj SK, Kahook MY, Gallardo MJ et al. Twelve-month outcomes of stand-alone excisional goniotomy in mild to severe glaucoma. Clin Ophthalmol. 2020; 14:1891-97.

ElMallah et al conducted a retrospective, multicenter case series of the KDB Glide® performed as a standalone procedure. Data was collected from 8 surgeons at 8 study sites (6 in the U.S. and 2 in Mexico). The purpose of the study was to characterize the reductions in intraocular pressure (IOP) and IOP-lowering medication following excisional goniotomy with KDB Glide® standalone in eyes with medically uncontrolled glaucoma. The study enrolled 42 eyes of 35 subjects with no restrictions for glaucoma sub-type or severity. The majority of eyes had a diagnosis of primary open angle glaucoma (n=36, 85.7%), but the study also included a small number of eyes diagnosed with pigmentary glaucoma (n=3, 7.1%), pseudoexfoliative glaucoma (n=2, 4.8%), and angle-closure glaucoma (n=1, 2.4%). Glaucoma severity was evaluated according to the International Classification of Diseases-10 (ICD-10) criteria, and was characterized as mild (19.0%), moderate (61.9%), or severe (19.0%). Subjects were followed for 12 months and were evaluated 1 day, 1 week, and 1, 3, 6, and 12 months postoperatively.

Baseline characteristics were as follows:

Baseline characteristics	KDB-Standalone (n=42)
Mean Age $\pm$ SE	71.3 $\pm$ 1.8 years
Sex	Female: 25 (59.5%) Male: 17 (40.5%)
Ethnicity	White: 19 (45.2%) Black: 13 (31.0%) Hispanic: 9 (21.4%) Asian: 1 (2.4%)
Medicated IOP Mean $\pm$ SE	21.6 $\pm$ 0.36 mmHG
Glaucoma Medications Mean $\pm$ SE	2.55 $\pm$ 0.22 meds

Included subjects were ≥18 years of age, phakic or pseudophakic, and with any glaucoma sub-type or severity. There were no baseline IOP requirements other than inadequate IOP control on at least 1 and up to 3 ocular hypotensive medications. Indication(s) for surgery included reducing IOP, reducing medication burden, or both. Subjects with a recent history (within 3 months) of IOP-lowering medication addition, laser trabeculoplasty, iridotomy, or initiation of systemic beta-blocker therapy were excluded from the study. Other exclusion criteria included any previous glaucoma surgery, and uncontrolled systemic conditions that might confound study measurements. Both eyes of a single subject were included if both qualified.

Following entry into the anterior chamber through a peripheral corneal incision, the KDB Glide® distal tip was used to pierce the trabecular meshwork (TM) and enter Schlemm's canal. The tip was advanced along the canal to elevate and guide the TM onto two parallel blades which then excised a strip of tissue. Approximately 3-5 clock hours of TM were excised.

The primary effectiveness endpoint was IOP reduction compared to baseline. The reduction in IOP-lowering medications was evaluated as the secondary effectiveness endpoint. Eyes with glaucoma secondary surgical intervention (SSI) for IOP control were excluded from the primary and secondary effectiveness analyses. Of the 42 eyes treated with KDB as a standalone procedure, 35 were available for the 12-month visit. Six eyes (n=6) required SSI for IOP control through 12 months and one eye (n=1) was considered lost to follow-up. At the 12-month visit, mean (SE) IOP was reduced by 3.94 (1.1) mmHg. Mean (SE) number of IOP-lowering medications was reduced by 0.31 (0.23) at 12 months (n=35).

In addition to the primary and secondary endpoints, ElMallah et al report additional effectiveness outcomes for all treated eyes (n=42) where eyes requiring SSI (n=6) and the eye that was lost to follow-up (n=1) were counted as failures. These additional analyses included the proportion of eyes with an IOP reduction of 20% or greater compared to baseline which was 42.9% (18/42 eyes).

Safety reporting consisted of the incidence of intraoperative and postoperative adverse events. The following adverse events were reported:

Adverse Event	KDB-Standalone (n=42)
Increased IOP	3 eyes (7.1%)
Descemet's Tear (intraoperative)	1 eye (2.4%)

The most common adverse event was postoperative elevated IOP (n=3, 7.1%). One case resolved with medical therapy, and the other two required a secondary surgical intervention. One eye experienced an intraoperative tear in Descemet's membrane, which was localized and did not require subsequent treatment during the 12-month study. No other adverse events were reported. Over the course of 12 months, 6 eyes required additional surgeries for inadequate IOP control. These included trabeculectomy (n=4; 2 of which had concurrent Ex-Press mini-shunt) and Ahmed valve implantation (n=2).

8.

Figure a

Figure b





Figure c

#### USE OF SYMBOLS

The following symbols used in the KDB GLIDE® packaging are defined below.

## 10.

	<b>English (en)</b>
	Manufacturer
	Date of Manufacture YYYY-MM
	Distributor
	"Use-by" or Expiration date
<b>QTY</b>	Quantity of devices in package.
<b>LOT</b>	Lot Number
<b>REF</b>	Model/Catalogue Number
<b>SN</b>	Serial Number
<b>MD</b>	Medical Device KDB GLIDE® Ophthalmic Instrument
	Single Sterile Barrier System
	Single Sterile Barrier System with Protective Packaging Outside
<b>STERILE R</b>	Sterilized using irradiation
	Do not use if package is damaged
	Do not re-sterilize
	Do not re-use
	Consult Instructions for Use or Consult Electronic Instructions for Use <a href="https://www.newworldmedical.com/product-ifus/">https://www.newworldmedical.com/product-ifus/</a>
	Caution
<b>Rx only</b>	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician
<b>UDI</b>	Unique Device Identifier
	Peel Direction
	Keep Dry
	Keep away from sunlight
<b>EC REP</b>	Authorized representative in the European Community

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Emergo Europe B.V.