



Ancillary Products

Patch Grafts

Pericardium: Freeze-Dried or Hydrated

Features:

- Biocompatible human tissue to cover the drainage implant tube
- Gamma sterilized
- Shelf life of 2 years
- Nominal thickness of 500 microns
- Size: 1.0 cm x 1.0 cm
- Available freeze-dried or hydrated
- Can be used as ocular tissue for other cosmetic uses
- Storage conditions:
 - Freeze-dried: Store at ambient temperature or cooler
 - Hydrated: Store at ambient temperature down to refrigerated temperature. Do not freeze



Sclera

Features:

- Biocompatible human tissue to cover the drainage implant tube
- Gamma sterilized
- Shelf life of 2 years
- Nominal thickness of 500 microns
- Size: 0.7 cm x 0.7 cm
- Can be used as ocular tissue for other cosmetic uses
- Storage conditions: Store between 2°C and 30°C. Do not freeze
- No rehydration required
- No rinsing needed




Tube Inserter™



Features:

- Notched tip secures the valve tube
- Provides rigidity to the valve tube for easy insertion into the anterior chamber
- Stainless steel
- Serrated grip
- Shorter tip helps to provide a better view through a microscope
- Reusable and inexpensive
- Can be used with any drainage device

Tube Extender™



Features:

- Provides extra tube length
- Provides a leak-proof junction between tubes
- Helpful when drainage tube is cut too short
- Can be used with other drainage devices
- Easily sutured onto the sclera

Ordering Information

TSH	Human Sclera, Hydrated
TPH	Pericardium, Hydrated
TPD	Pericardium, Freeze-Dried
Model: TI	Tube Inserter
Model: TE	Tube Extender



New World Medical
10763 Edison Court
Rancho Cucamonga, CA 91730
(909) 466-4304

www.newworldmedical.com

For more information or to place an order, contact your local representative.